BioDtex™

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 biuro@interarma.pl


 www.interarma.pl/BioDtex

BioDtex

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BioDtex[™]

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Report on:
Assessment of BioDt看ex unit instrument for detection of food industry biofilms
Work performed by Campden BRI (Chipping Campden) Limited
Report number: MB/REP/ZCL/156947/1 ￼ Issue date: 26th September 2022

Contact details:
Zoe Lambert ￼ Microbiology ￼ Campden BRI (Chipping Campden) Limited
Zoe.Lambert@campdenbri.co.uk ￼ Tel: +44(0)1386 842274

Report issued and authorised by:
Campden BRI (Chipping Campden) Limited
Rob Limburn ￼ Section Lead, Industrial Process Microbiology

Our ref: MB/REP/Z/156947/1
Page count: 19

￼ ￼ ￼ ￼

￼ Prior written approval is required for use of the Campden BRI logo and for the reproduction of this report if not reproduced in full.

￼ Unless this report includes an explicit statement of compliance/non-compliance with a requirement and/or specification, no such statement should be inferred.

￼ Unless this report includes an explicit statement to the contrary, results reported relate only to the items tested.

￼ The information provided within this document is given after the exercise of all reasonable care and skill in its compilation, preparation and issue but is provided without liability in its application and use. Any opinions and interpretations are not provided under the auspices of any third-party certification or accreditation. Unless otherwise expressly agreed in writing and signed by a duly authorised representative of Campden BRI (Chipping Campden) Limited, this services was subject to our Standard Terms and Conditions of Contract, available on request or from our website: <http://www.campdenbri.co.uk/campdenbri/CampdenBRISupplyTerms.pdf>

The information in this document is only intended for the individual or entity to whom it is addressed. It may contain privileged and confidential information that is exempt from disclosure by law and if you are not the intended recipient, you must not copy, distribute or take any action in reliance on it. If you have received this document in error please notify us immediately by telephone on +44(0)1386 842000.

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